

Powered by



Powered by



Hoof Link System

QH Chip

Horse is implanted with an ISO approved low frequency microchip. The HoofLink System recommends thermal microchips placed in the lip, but allows other placements.

Authorized users scan the horse and use the microchip ID to enroll their details in the HoofLink Application.

AQHA identifier scans the horse, and enters the microchip number into the QH Chip Application.

Enrollment includes ownership information, description of markings, health documents, and four ID photos. Data is immediately uploaded to the HoofLink Database.

AQHA Identifier inspects and takes four photos of the horse and uploads the photos into the QH Chip Application.

Payment is taken by cash, check, or credit card.

The transaction is finalized through payment by credit card.

Horse is now enrolled, horsemen can now scan their horses and take daily temperatures, view health documents, and check in and out of equine venues.

Horse is now officially verified. The AQHA identifier stamps the horse's registration certificate and applies the microchip barcode, if available.

Horsemen can view their data anytime from the secure database. Additional users can view this data only with the permission of the horsemen.

The official identification details are transmitted to AQHA's database immediately, which feeds to Incompass Solutions within 24 hours.

The HoofLink System is used in the paddock, test barn, stable gate, and during works to identify horses and more.

The QH Chip Application is intended for internal use only by AQHA staff and identifiers for verifying horses prior to racing.

The HoofLink System is paired with advanced Bluetooth Scanners provided by Lip Chip LLC.

Lip Chip is a corporate partner with the AQHA. The QH Chip application utilizes microchip scanners provided by Lip Chip built specifically for the AQHA.